

**Stratus OCT™** with Software Version 5.0

**Real Answers in Real Time**



# Stratus OCT™ — The standard of care system for comprehensive retina and glaucoma management.




**Visualize virtual histology.** High resolution cross-sectional imaging for the comprehensive management of glaucoma and retinal disease.

**Assess RNFL change.** GPA Advanced Serial Analysis reports statistically significant change and rate of change in RNFL thickness.

**Expand your clinical confidence.** Most extensive library of clinical studies in the industry and over 8000 systems installed worldwide.

**Offer comprehensive care.** Valuable for pre- and post-op cataract patients to identify and illustrate cause of poor vision.


**Increase practice productivity.** Operate the full analysis capabilities from the lane, your office or a remote location.




In the Stratus OCT image display, retinal layers with the highest reflectivity appear red. In a healthy retina, these include the nerve fiber layer, retinal pigment epithelium and choriocapillaris. The layers that exhibit minimal reflectivity appear blue or black, such as the photoreceptor layer, choroid, vitreous fluid or blood.

# Obtain real-time non-invasive histology of live tissue


Stratus OCT reveals the retinal layers in high-resolution, cross-sectional views, offering insight for diagnosis, therapy and ongoing management of retinal disorders.




Cystoid Macular Edema




Central Serous Chorioretinopathy




Age-related Macular Degeneration with Overlying Cystoid Macular Edema




Macular Hole with Operculum




Diabetic Macular Edema



Epiretinal Membrane with Lamellar Hole and Cystoid Macular Edema




Branch Retinal Vein Occlusion




RNFL Loss

# Visualize and analyze retinal disorders




## Diabetic Retinopathy

- Stratus OCT reveals and measures diffuse macular thickening and loss of foveal contour
- Intraretinal cysts and fluid accumulation are identifiable as areas of low reflectivity in the cross-sectional scan
- Post-treatment resolution of retinal thickening can be quantified and monitored



## Epiretinal Membrane


- Stratus OCT scan shows the epiretinal membrane as a highly reflective band on the inner retinal surface
- Separation of the membrane from the retina is visible in areas
- Underlying retina is thickened, with loss of normal foveal contour



## Age-related Macular Degeneration


- Disruption of RPE, caused by neovascularization and drusen, can be visualized
- Pockets of interretinal fluid are visible as areas of reduced reflectivity
- Structural changes resulting from therapy can be quantified and monitored

# Detect glaucoma damage at an earlier stage




## Retinal Nerve Fiber Layer Analysis

- Analysis of RNFL aids in identification of early glaucomatous loss
- Circular scans of 3.4 mm diameter around optic nerve head provide measurement of RNFL in the peripapillary region
- RNFL thickness measurement is graphed in a TSNIT orientation and compared to age-matched normative data



## Optic Nerve Head Analysis

- Radial line scans through optic disc provide cross-sectional information on cupping and neuroretinal rim area
- Disc margins are objectively identified using signal from end of RPE
- Key parameters include cup-to-disc ratio and horizontal integrated rim volume<sup>1</sup>



## Macular Thickness Analysis

- Thinning of the macula may reflect glaucomatous loss
- Structural analysis of retinal sublayers reveals macular complications
- Cross-sectional view provides visualization and measurement of retinal layers

# Stratus OCT Printout

## Retinal Thickness Report


**Scan Protocol:** Fast Macular Thickness, Macular Thickness, Line, Cross Hair  
**Used for:** Assessment of overall macular region or specific areas of interest

**STRATUS OCT**  
Retinal Thickness Report - 4.0.3 (0070)


DOB: 7/9/1924, ID: NA, Gender

Scan Type: Fast Macular Thickness Map: OS  
 Scan Date: 3/31/2004  
 Scan Length: 6.0 mm

**OCT Image**



**Fundus Image**



Signal Strength (Max 10) 6

Retinal Thickness is 451 microns at A-scan 55

Caliper Length is OFF


**Legend**

Arrow indicates direction and angle of individual scan displayed in scan image.

**Normative Data**

Distribution is indicated by a stoplight color code, described in detail inside back cover.

**Thickness Chart**



Signature: \_\_\_\_\_

Physician: \_\_\_\_\_

**Patient Information** →

**Scan Information** ←

**Scan Image** →  
The image with the lowest signal strength or with a data message will be shown to facilitate operator review.

**Fundus Image** ←  
Scan should be centered on fovea.

**Signal Strength** →  
Numbers range from 0 (weak) to 10 (strong). Analysis algorithm may fail on scans with low signal strength; therefore, images with a signal strength below 5 should not be used for analysis.

**Legend** →

**Normative Data** ←

**Data Message** →  
"Scan Too High," "Scan Too Low" or "Missing Data" message will be shown, if applicable, regarding placement and completeness of scan.

**Thickness Chart** →  
Graphic display of retinal thickness. Colored bands demonstrate range of normative data.

**Retinal Thickness** →  
Displayed for A scan indicated. This corresponds to A scan (shown here) selected on thickness chart.

**Caliper Length** →  
Measurement indicates distance between calipers if they are placed on scan image during analysis.

# Retinal Thickness Tabular Output

Scan Protocol: Radial Lines, Fast Macular Thickness, Macular Thickness

Used for: Imaging and measurement of macular pathology

**Patient Information**

DOB: 7/9/1924, ID: NA, Gender:

**Scan Image**

The image with the lowest signal strength or with a data message will be shown to facilitate operator review.

**Fundus Image**

Scan should be centered on the fovea.

**Signal Strength**

Numbers range from 0 (weak) to 10 (strong). Analysis algorithm may fail for scans with low signal strength; therefore, images with a signal strength below 5 should not be used for analysis.

**Data Message**

"Scan Too High," "Scan Too Low" or "Missing Data" message will be shown, if applicable, regarding placement and completeness of scan.

**Map**

Thickness is displayed using a color scale. A normal eye will be displayed as blue centrally. Orange and red indicate greater thickness. Legend for color scale appears near bottom right of page.

**Sector Averages**

Numerical values indicate average thickness of each sector. Colors within each sector indicate comparison to normative data.

**Scans Used**

If any of the 6 radial line scans has been deselected, analysis might be incomplete.

**STRATUS OCT**  
Retinal Thickness Tabular Output Report - 4.0.3 (0070)

Scan Type: Fast Macular Thickness Map  
Scan Date: 3/31/2004  
Scan Length: 6.0 mm

DOB: 7/9/1924, ID: NA, Gender:

OD OS

Map Diameters: 10mm, 30mm, 60mm

Parameter	OD	OS	Diff (OD-OS)
<b>Thickness</b>			
Foveal minimum	129	429	-290
Fovea	169	425	-256
Temporal inner macula	226	327	-101
Superior inner macula	226	384	-148
Nasal inner macula	236	354	-118
Inferior inner macula	237	336	-99
Temporal outer macula	184	232	-48
Superior outer macula	184	251	-65
Nasal outer macula	220	204	+14
Inferior outer macula	193	220	-27
Superior/Inferior outer	1.016	1.091	-0.075
Temporal/Nasal inner	0.954	0.850	0.102
Temporal/Nasal outer	0.836	0.870	-0.043
<b>Volume (cubic mm)</b>			
Fovea	0.133	0.333	-0.200
Temporal inner macula	0.355	0.514	-0.159
Superior inner macula	0.372	0.603	-0.231
Nasal inner macula	0.371	0.62	-0.249
Inferior inner macula	0.373	0.527	-0.154
Temporal outer macula	0.375	1.204	-0.259
Superior outer macula	1.042	1.322	-0.260
Nasal outer macula	1.171	1.413	-0.232
Inferior outer macula	1.025	1.22	-0.195
Total macula volume	5.020	7.790	-1.969

Parameter Normal Range

Fovea Minimum: Measurement at center of fovea where radial scan lines intersect	135 – 215 µm
Average Thickness	
Fovea	168 – 239 µm
Temporal Inner Macula	240 – 294 µm
Superior Inner Macula	243 – 296 µm
Nasal Inner Macula	240 – 297 µm
Inferior Inner Macula	246 – 297 µm
Temporal Outer Macula	199 – 276 µm
Superior Outer Macula	207 – 256 µm
Nasal Outer Macula	198 – 274 µm
Inferior Outer Macula	207 – 256 µm
Ratios	
Superior/Inferior Outer	0.832 – 1.222
Temporal/Nasal Inner	0.800 – 1.227
Temporal/Nasal Outer	0.557 – 1.845
Volume	
Fovea	0.13 – 0.19 cubic mm
Temporal Inner Macula	0.38 – 0.46 cubic mm
Superior Inner Macula	0.38 – 0.46 cubic mm
Nasal Inner Macula	0.38 – 0.47 cubic mm
Inferior Inner Macula	0.39 – 0.47 cubic mm
Temporal Outer Macula	1.06 – 1.46 cubic mm
Superior Outer Macula	1.10 – 1.36 cubic mm
Nasal Outer Macula	1.05 – 1.45 cubic mm
Inferior Outer Macula	1.10 – 1.35 cubic mm
Total Macula Volume	6.18 – 7.42 cubic mm

Normal distribution, macula thickness normative data, Carl Zeiss Meditec.

**Legend**

Normative data is displayed in spotlight color code, described in detail inside back page.

**Color Scale**

For thickness maps.

# Stratus OCT Printout

## Optic Nerve Head Analysis Report

Scan Protocol: Optic Disc, Fast Optical Disc

Used for: Evaluation of the optic disc

### Patient Information

#### Cup Markers

The edge of the cup is indicated with a green dot on the scan image and a green x on the composite diagram.

#### RPE Markers

The end of the RPE is indicated. Shown blue on the scan image, red on the composite diagram.

#### Signal Strength

Numbers range from 0 (weak) to 10 (strong). Analysis algorithm may fail on scans with low signal strength; therefore, images with a signal strength below 5 should not be used for analysis.

#### Optic Nerve Head Results

Data are derived from the 6 radial line scans.

#### Vertical Integrated Rim Area

Total volume of RNFL tissue in the rim is obtained by multiplying the average of the 6 individual rim areas by the circumference of the disc. Normal values are  $0.36 \pm 0.08$  cubic mm.<sup>2</sup>

#### Horizontal Integrated Rim Width

Total rim area is obtained by multiplying the average of the 6 individual nerve widths by the circumference of the disc.

#### Disc Area

The area within the red outline on the composite diagram.

#### Cup Area

The area within the green outline on the composite diagram.

#### Rim Area

Disc area minus cup area.

#### Cup/Disc Area Ratio

Ratio of cup area to disc area.

#### Cup/Disc Horizontal Ratio

Ratio of the longest horizontal line across the cup to the longest horizontal line across the disc.

#### Cup/Disc Vertical Ratio

Ratio of the longest vertical line across the cup to the longest vertical line across the disc.

#### Fundus Image

Can be used to verify scan placement.

#### Composite Diagram

Yellow line indicates individual scan selected and displayed above.

#### Scan List

For each of the 6 radial line scans, this indicates the date when the scan was last adjusted and saved.

#### Cup Area


The area within the green outline on the composite diagram.

#### Cup Volume

Volume of cup area on the topographic display.

#### Cup Offset for Topo

Default is 150  $\mu$ m anterior to disc. Represents edge of blue shading on topographic display—the surface of the cup as if it were filled with water.



### Scan Information

#### Individual Radial Scan Analysis

Data relates only to the individual scan image displayed.

#### Rim Area

Indicated with red shading on scan image, this area is bounded by the cup diameter line and a line from the RPE marker to the anterior surface of the disc, at a 90-degree angle to the cup diameter line.

#### Average Nerve Width @ Disc

The average of the nerve bundle widths at the disc on each side. Nerve bundle width indicated by yellow line from RPE marker to anterior surface.

#### Disc Diameter

Illustrated and measured on a straight line between the 2 RPE markers. Blue line on scan image.

#### Cup Diameter

Illustrated and measured on a straight line parallel to, and 150  $\mu$ m anterior to, the disc diameter line. Red on image, green on diagram, this line is adjustable.

#### Horizontal Rim Length

Disc diameter minus the cup diameter.

#### Legend

Arrow indicates direction and angle of individual scan displayed in scan image.

#### Cup Offset

Point at which cup diameter is measured. Default placement is 150  $\mu$ m anterior to the RPE markers; this point is adjustable.

#### Plot Background

There are 3 optional backgrounds; two are similar to a topographic map. **Aligned and Shaded** (default) provides a shaded relief topographic map and corrects for patient motion. **Absolute** provides a non-shaded topographic map with no correction for movement. **None** shows a solid gray background with no appearance of depth.



# RNFL Thickness Average Analysis

Scan Protocol: RNFL 3.4 mm, Fast RNFL 3.4 mm

Used for: Retinal nerve fiber layer thickness assessment and comparison to normative database

**Patient Information**

DOB: 9/9/1953, ID: NA, Female

**RNFL Thickness Chart**  
Peripapillary RNFL thickness is displayed in TS-NIT format. Colored bands demonstrate range of normative data.

**Sector Averages**  
Comparison to normative data in each sector is indicated with stoplight color scheme. Values are displayed numerically.

**Quadrant Averages**  
Comparison to normative data in each quadrant is indicated with stoplight color scheme. Values are displayed numerically.

**OD/OS Graph**  
TS-NIT line graph displays RNFL thickness in both eyes. Asymmetry may be indicative of glaucomatous loss.

**Scans Used**  
If any of the 3 circular scans has been deselected, analysis might be incomplete and for the fast RNFL scan, normative data will not be displayed.

**Legend**  
Normative data is displayed in stoplight color code, described in detail on inside back cover.

**Scan Information**

Scan Type: Fast RNFL Thickness (3.4)  
Scan Date: 10/29/2003  
Scan Length: 10.87 mm

**Signal Strength**  
Numbers range from 0 (weak) to 10 (strong). Analysis algorithm may fail on scans with low signal strength; therefore, images with a signal strength below 5 should not be used for analysis.

**Data Message**  
"Scan Too High," "Scan Too Low" or "Missing Data" message will be shown, if applicable, regarding placement and completeness of scan.

	OD (N=3)	OS (N=3)	OD-OS
lmax/Smax	1.09	1.21	-0.12
Smax/lmax	0.92	0.83	0.09
Smax/Tavg	2.28	1.65	0.73
lmax/Tavg	2.49	1.68	0.81
Smax/Navg	1.78	2.08	-0.30
Max-Min	150.00	96.00	54.00
Smax	179.00	105.00	74.00
lmax	195.00	127.00	68.00
Savg	147.00	71.00	76.00
lavg	150.00	76.00	74.00
Avg.Thick	118.66	66.22	52.45

**Tabular Data**  
For each eye, values are compared to normative data and displayed in stoplight color code, described in detail on inside back cover.

Comparison	Average Ratios	Measurement	Average Values
<b>lmax/Smax</b> Thickest points in inferior and superior quadrants	0.80 – 1.25	<b>Min-Max</b> Difference between minimum and maximum measurements	96 – 154 µm
<b>Smax/lmax</b> Thickest points in superior and inferior quadrants	0.77 – 1.25	<b>Smax</b> Thickest measurement in superior quadrant	124 – 189 µm
<b>Smax/Tavg</b> Thickest point in superior quadrant to average in temporal	1.70 – 3.06	<b>lmax</b> Thickest measurement in inferior quadrant	125 – 194 µm
<b>lmax/Tavg</b> Thickest point in inferior quadrant to average in temporal	1.69 – 3.12	<b>Savg</b> Average measurement in superior quadrant	97 – 152 µm
<b>Smax/Navg</b> Thickest point in superior quadrant to average in nasal	1.37 – 2.93	<b>lavg</b> Average thickness in inferior quadrant	98 – 156 µm
		<b>Average Thickness</b>	82 – 118 µm

Normal distribution, RNFL normative data, Carl Zeiss Meditec.

# Stratus OCT Printout

## GPA™ Advanced Serial Analysis

Scan Protocol: Fast RNFL Thickness (3.4), RNFL Thickness (2.27 x disc)

Used for: Statistical analysis of RNFL thickness change over time. Can be applied to up to 8 OD and/or 8 OS scan groups

**Patient Information**

Dx, GPA Change OD No Change OS

DOB: 1/1/1960, ID: GPA0007, Female

**Scan Information**

Scan Type: Fast RNFL Thickness (3.4)

Scan Date: Multiple, See Table

Scan Length: 10.87 mm

**TSNIT Graph**

Peripapillary RNFL thickness is displayed in TSNIT format. Line colors correspond to scan date, as shown in the RNFL summary table.

**RNFL summary table**

	OD	SS,Q	AVG	SUP	INF
3/13/2003 11:16:07 AM (N=2)	8	87.80	112.00	100.00	
5/20/2004 10:46:57 AM (N=2)	8	75.24	103.00	112.00	
6/15/2005 10:19:55 AM (N=2)	8	71.25	104.00	111.00	
1/6/2006 10:24:25 AM (N=2)	7	84.52	87.00	116.00	
1/4/2007 6:55:39 AM (N=2)	8	80.46	95.00	114.00	

Rate of change:  $-5.926 \pm 1.194 \mu\text{year}^*$   
Statistically significant  $P < 0.1\%$ , seek clinical correlates

**Thickness over time graph**

Average RNFL thickness for each exam is plotted. Dot colors correspond to scan date, as shown in the RNFL summary table.

**Rate of change**

	OS	SS,Q	AVG	SUP	INF
3/13/2003 11:16:07 AM (N=2)	8	87.80	116.00	124.00	
5/20/2004 10:46:27 AM (N=2)	7	87.80	104.00	116.00	
6/15/2005 10:20:53 AM (N=2)	8	96.83	114.00	124.00	
1/6/2006 10:22:29 AM (N=2)	8	86.83	116.00	104.00	
1/4/2007 9:59:42 AM (N=2)	8	78.30	112.00	88.00	

Rate of change:  $-3.189 \pm 7.58 \mu\text{year}^*$   
Statistically not significant  $P > 5\%$

**Thickness over time graph**

**Legends**

\* 95% confidence interval

QC	Normal Distribution Percentile
A - Scan too low	100%
B - Scan too high	95%
C - Low confidence in analysis	5%
D - Missing data	1%
E - Edited Layer	0%
F - No Normative Data	0%

Signature: \_\_\_\_\_


Physician: \_\_\_\_\_

Site ID: ZEISS Clinic

# Stratus OCT Normative Data


## Stoplight Color Scheme

RNFL Normative Data Display




RNFL Normative Distribution

Of the normal population:




Macula Normative Data Display



Macula Normative Distribution

Of the normal population:



At Carl Zeiss Meditec, we are committed to expanding clinical potential with innovative, precise and clinically advanced instruments that contribute to the enhancement of vision worldwide. And, through lifetime customer care and ongoing technology upgrades, we are also dedicated to ensuring your success now and throughout the future. For more information on the Stratus OCT system or to order, contact your Carl Zeiss Meditec representative today, or visit our website at [www.meditec.zeiss.com/stratus](http://www.meditec.zeiss.com/stratus).

1. Wollstein G, Ishikawa H, Wang J, Beaton SA, Schuman JS. Comparison of three optical coherence tomography scanning areas for detection of glaucomatous damage. *Am J Ophthalmol*. 2005;139(1):39-43.
2. Schuman JS, Wollstein G, Farra T, et al. Comparison of optic nerve head measurements obtained by optical coherence tomography and confocal scanning laser ophthalmoscopy. *Am J Ophthalmol*. 2003;135(4):504-512.

Pentium is a registered trademark of Intel Corp.  
Windows is a registered trademark of Microsoft Corp.

## Technical Specifications

Tomographic Imaging	
Purpose	Cross-sectional imaging of retina
Signal type	Optical scattering from tissue
Signal source	Superluminescent diode, 820 nm
Optical power	≤750 microwatts at cornea
Longitudinal/Axial resolution	≤10 µm in tissue
Transverse sample size	20 µm in tissue
Scanners	Galvanometric mirror
Scan patterns	Line, circle, concentric rings, radial lines
Scan pixels	Adjustable from (1024 axial x 128 transverse) to (1024 axial x 768 transverse)
Longitudinal (depth) range	2 mm in tissue
Scan rate	400 A scan/sec
Fundus Imaging	
Purpose	Fundus alignment, documentation
Signal type	CCD image
Field of view	26° x 20.5°
Viewing method	Flat panel display
Illumination	Near IR/red-free
Internal fixation	32 x 16 LED dot matrix
External fixation	Slit lamp type adjustable blinking LED
Minimum pupil diameter	3.2 mm
Electrical	
Power consumption	100 V approx. (±10%), 50/60 Hz, 6.0 A 115 V approx. (±10%), 60 Hz, 6.0 A 230 V approx. (±10%), 50/60 Hz, 3.0 A 700 VA
Footprint	
Patient module	48 inches x 34 inches, 120 cm x 85 cm
User Features	
Processor	2.4 GHz Pentium® IV
Operating system	Windows® 2000
Memory	512 MB
Standards and Approvals	
UL 2601-1	
CSA 22.2 No. 601.1	
MDD	

Note: All technical specifications are subject to change without notice.

### Carl Zeiss Meditec AG

Goeschwitzer Str. 51-52  
07745 Jena  
Germany  
Telefon: +49 (0)36 41 / 2 20-3 33  
Telefax: +49 (0)36 41 / 2 20-2 82  
[info@meditec.zeiss.com](mailto:info@meditec.zeiss.com)  
[www.meditec.zeiss.com](http://www.meditec.zeiss.com)

### Carl Zeiss Meditec, Inc.

5160 Hacienda Drive  
Dublin, CA 94568  
USA  
Tel: 1-925-557-4100  
Toll free: 1-800-342-9821  
Fax: 1-925-557-4101  
[info@meditec.zeiss.com](mailto:info@meditec.zeiss.com)  
[www.meditec.zeiss.com](http://www.meditec.zeiss.com)